

**SOP#112: Agencourt DNA extracdction from saliva colleted in Genotek Oragene kits**

**1. PURPOSE:**

This protocol describes the manual purification of DNA using Agencourt’s DNAdvance kit for saliva collected using Oragene® self-collection kits.

**2. MATERIALS**

**2.1 Reagents and supplies**

<b>Material</b>	<b>Vendor</b>	<b>Cat#</b>
Saliva samples collected from Oragene self-collection kit		
1.5mL microcentrifuge tubes	Varies	Varies
1.2mL storage plate	Thermo Fisher	AB-1127
P1000 pipetteman and filter tips	Varies	Varies
P200 pipetteman and filter tips	Varies	varies
Agentcourt SPRIPlate Super magnet	Alpaqua	A32782
Bind1 buffer	Beckman-Coulter Genomics	Part of Agencourt DNAdvance kit (Cat# A48705)
Bind2 buffer	Beckman-Coulter Genomics	Part of Agencourt DNAdvance kit (Cat# A48705)
Elution buffer	Beckman-Coulter Genomics	Part of Agencourt DNAdvance kit (Cat# A48705)
70% Ethanol	Varies	Varies
Hard-Shell® 96-Well PCR Plates	Bio-Rad	HSP9901
Agarose	Varies	Varies

**2.2 Equipment**

- 50°C water bath
- 4°C fridge or cold room
- 80°C freezer
- Nanodrop 2000 Spectrometer
- Gel boxes
- Power supply
- Gel imager

### 3. EXPERIMENTAL PROCEDURES

#### 3.1 Sample Preparation

- 3.1.1. Collect saliva samples with the Oragene self-collection kit. Samples can be stored at room temperature for at least 3 years.
- 3.1.2. Incubate samples at 50°C for 1 hour. Samples can be stored at room temperature for at least 3 years after this step.
- 3.1.3. Enter the sample IDs to the 96 well plate map and print a hard copy.
- 3.1.4. Transfer 500 µL of saliva sample into a 1.2mL storage plate according to the plate map.

#### 3.2 Binding

- 3.2.1. Add 200 µL of Bind1 buffer to each well and mix by pipetting 10 times or until mixed well.
- 3.2.2. Shake Bind2 buffer bottle to until the bead particles are re-suspended well.
- 3.2.3. Add 340 µL of Bind2 buffer in each well and mix by pipetting 15 times or until mixed well. Pipette slowly to avoid air bubbles.
- 3.2.4. Incubate the plate at room temperature for 1 minute.
- 3.2.5. Place the sample plates on an Agencourt SPRIPlate Super magnet for 8 minutes to separate.
- 3.2.6. Aspirate and discard the supernatant while the plate is still situated on the magnet. Avoid disturbing the ring of magnetic beads; it's okay to leave the last 10 – 20 µL.

#### 3.3 Ethanol wash

- 3.3.1. Remove plate off magnet. Add 700 µL of 70% ethanol to each well and mix by pipetting 20 times to re-suspend beads from the bottom of the well. Try to eliminate any bead clumps with tip mixing.  
**NOTE:** Make fresh 70% ethanol for each extraction.
- 3.3.2. Place the plate back on the magnet for 2 minutes, or until the solution clears.
- 3.3.3. Aspirate and discard the supernatant while the plate is still situated on the magnet. Avoid disturbing ring of magnetic beads.
- 3.3.4. Repeat steps 3.3.1 through 3.3.3 two more times.
- 3.3.5. Remove as much of the final ethanol wash as possible before adding the elution buffer.

#### 3.4 Elution

- 3.4.1. Remove the plate off the magnet. Add 50 µL of elution buffer to each well and mix by pipetting 10 times or until the magnetic beads are re-suspended completely from the bottom of the well.
- 3.4.2. Place the plate back on the magnet for 3 minutes, or until the solution clears.

#### 3.5 Storage

- 3.5.1. Label a clean 96-well plate as the source DNA plate for the specific study. Transfer 40 µL of supernatant to the labeled 96-well plate for storage. Slowly aspirate to avoid disturbing the ring of beads; If beads are aspirated during transfer, dispense the sample back into the plate and incubate for another 5 minutes before repeating aspiration.
- 3.5.2. Quantify DNA by measuring OD260/OD280. Enter the DNA concentration into the sample manifest and calculate the amount of source DNA and elution buffer needed to make the normalized plate, based on the template.
- 3.5.3. Label a new clean 96-well plate as the normalized DNA plate. Normalize DNA to 40 ng/ul in a new 96-well plate, with well positions matching the source DNA plate.

**[Blackburn Lab SOP#112: Agencourt DNA extraction from saliva collected in Genotek Oragene kits Rev B]  
[Jan 2019]**

3.5.3. Check DNA integrity by running 3 ul of the normalized DNA on 0.8% agarose gels.  
Exclude degraded samples in the final telomere length analysis report.  
3.5.2. Store DNA at -80°C for further use.

**REVISION HISTORY**

DATE	AUTHOR	REVISION #	REVISION REASON
March 13, 2018	Calvin Wu	A	Initial Release
Jan 8, 2018	Jue Lin	B	Added DNA quantification and quality check. A few other minor revisions

FINAL JAN 2019